

Dental Chair - Check-list Before You Order

1. Is the dental chair's **floor control box** to be positioned inside (built-in) or outside (external to) the chair? [Floor Control Box / Check-list Before You Order / Dental Chairs - Documents & Information / Documents & Information | Al-Mahdi - Dental Supplies](#)
2. Does your surgery have appropriate suction/vacuum and compressor system for the chair? (**Amalgam separator** is required for UK suction systems). [Suction with amalgam separator / Check-list Before You Order / Dental Chairs - Documents & Information / Documents & Information | Al-Mahdi - Dental Supplies](#)
3. Do you have appropriate plumbing and electronics installed for the chair?
4. Is your surgery large enough for the **chair dimensions**? [Chair Dimensions / Check-list Before You Order / Dental Chairs - Documents & Information / Documents & Information | Al-Mahdi - Dental Supplies](#)
5. Is the chair intended for left handed or right handed use?
6. Do you need a UK power plug or is the chair intended for a different region?
7. Do you need 2 hole, 4 hole or 6 hole tubing for the hand piece attachments? (The chairs come with 4 hole tubing as standard)
8. Do you need fibre optic tubing for the highspeed attachment?

9. Do you need new high speed/slow speed handpieces, ultrasonic scalers and matching scaler tips?
10. Do you need any dental chair **add-ons**? [Any \(Dental Chair\) Add-ons / Dental Chairs / WEBSHOP | Al-Mahdi - Dental Supplies](#)
11. What's the name and contact details of your local dental equipment engineer (for servicing dental chairs).

Please send us an email with the answers to the check-list **PRIOR to placing an order.**

We want to make sure that you are happy with your order and that you receive a product that is beneficial for your practice and meets your clinical requirements.

Al-Mahdi Dental Supplies
www.al-mahdi-dental-supplies.com
Email: info@al-mahdi-dental-supplies.com
Tel & WhatsApp: +447770447965