

Dental Chair - Check-list Before You Order

1. Is the dental chair's **floor control box** to be positioned inside (built-in) or outside (external to) the chair? [Floor Control Box / Check-list Before You Order / Dental Chairs - Documents & Information / Documents & Information | Al-Mahdi - Dental Supplies](#)
2. Does your surgery have appropriate suction/vacuum and compressor system for the chair? ***Amalgam separator** is required for UK suction systems. Regulations in the UK requires that the chair's spittoon and suction water is undergoing amalgam separation. [Suction with amalgam separator / Check-list Before You Order / Dental Chairs - Documents & Information / Documents & Information | Al-Mahdi - Dental Supplies](#)
3. Do you have appropriate plumbing and electronics installed for the chair?
4. Is your surgery large enough for the **chair dimensions**? [Chair Dimensions / Check-list Before You Order / Dental Chairs - Documents & Information / Documents & Information | Al-Mahdi - Dental Supplies](#)
5. Is the chair intended for left handed or right handed use?
6. Do you need a UK power plug or is the chair intended for a different region?
7. Do you need 2 hole, 4 hole or 6 hole tubing for the hand piece attachments? (The chairs come with 4 hole tubing as standard)

8. Do you need fibre optic tubing for the highspeed attachment?

9. Do you need new high speed/slow speed handpieces, ultrasonic scalers and matching scaler tips?

10. What's the name and contact details of your local dental equipment engineer (for servicing dental chairs).

11. Have you read through, understood and accepted the terms and conditions stated in the "[Dental Chairs - Important info / Q & A](#)" page?

Please send us an email with the answers to the check-list **PRIOR to placing an order.**

We want to make sure that you are happy with your order and that you receive a product that is beneficial for your practice and meets your clinical requirements.

Al-Mahdi Dental Supplies
www.al-mahdi-dental-supplies.com
Email: info@al-mahdi-dental-supplies.com
Tel & WhatsApp: +447770447965