

Dental Chair - Check-list Before You Order

1. Does your surgery have an appropriate suction/vacuum and compressor system for the chair?
2. For UK customers: Regulations in the UK requires that all waste water (the chair's spittoon and suction water) is undergoing amalgam separation. Al-Mahdi Dental Supplies will ensure that the dental chair will be configured by the manufacturer so all the waste water tubing is collected in one pipe, which can then be connected to the suction system in the practice by your dental chair engineer. It's the customer's responsibility to ensure that the suction system in the practice is equipped with an [amalgam separator](#).
3. Do you have appropriate plumbing and electronics in the surgery prior to receiving the dental chair?
4. Is your surgery large enough for the chair dimensions?
5. Is the chair intended for left handed or right handed use?
6. Do you need a UK power plug or is the chair intended for a different region?
7. Do you need 2 hole, 4 hole or 6 hole tubing for the hand piece attachments? (The chairs come with 4 hole tubing as standard)
8. Do you need fibre optic tubing for the highspeed attachment?
9. Do you need new high speed/slow speed handpieces and ultrasonic scalers?
10. Can you confirm that you have a local dental equipment engineer who has agreed to install and service the chair if needed? (Al-Mahdi Dental Supplies will provide spare parts for free during the warranty period if needed)

Please send us an email with the answers to the check-list **PRIOR to placing an order.**

We want to make sure that you are happy with your order and that you receive a product that is beneficial for your practice and meets your clinical requirements.